

Rightest™ GS300 BLOOD GLUCOSE TEST STRIP INSERT

Intended Use

The **Rightest™** Blood Glucose Monitoring System is used by individuals with diabetes. It's for checking on glucose levels of whole blood from capillary. Capillary blood can be sampled from the fingertip, palm and forearm. It's as an aid in management of diabetes at home and clinical sites.


Rightest™ Blood Glucose Test Strips are intended for testing outside the body (*in vitro* diagnostic use) (For self-testing) only.

The **Rightest™** System tests the capillary blood and provides results equivalent to a laboratory instrument (plasma equivalent)

- The **Rightest™** Blood Glucose Test Strip GS300 is designed for use only with the **Rightest™** Blood Glucose Meter GM300.
- The **Rightest™** Blood Glucose Monitoring System includes Meter, Test Strips, Smart Code Key, Check Key, Control Solutions, Lancing Device and Lancets.

Test Procedure

REFER TO THE **Rightest™** USER'S MANUAL FOR MORE DETAILED INFORMATION.




Smart Code Key Installation


- 1) With the Meter off, follow ① and ② direction to put the new Smart Code Key into the track on code key base.
- 2) Push down the Smart Code Key until it snaps into the Smart Code Key base.

Preparing the Lancing Device


1. Pull off the depth adjustable cap.
2. Insert a new disposable lancet firmly into lancet carrier.
3. Twist off and set aside the protective cover of the disposable lancet.
4. Replace the depth adjustable cap.
5. Choose a depth of penetration by rotating the top portion of the depth adjustable cap until the setting depth matches the window. Settings are window. Settings are based on skin type "□□□" for soft or thin skin; "□□□□" for average skin; "□□□□□" for thick or calloused skin.
6. Hold the hub in one hand and pull on the plunger in the other hand. The device will be cocked. Release the plunger, it will automatically move back to its original position near the hub.




Performing a Test



- 1) Wash and dry your hands. Take one strip from the vial. Close the vial cap immediately.
- 2) Insert the strip into the strip port on meter with the indication symbol facing up. Push the strip in until it snaps and stops. The meter will turn on automatically.
- 3) Make sure the code number on the meter screen matches the code number on the test strip vial.
- 4) When you see the flashing blood drop, hold the lancet device to side of your fingertip and press the release button.
- 5) Gently squeeze your fingertip to get a drop of blood. Our meter only needs a tiny blood sample.





Sample Size Example



Please take a minimum of 1.4 µL to do the test on blood glucose monitoring system. Blood sample size above 4.0 µL might contaminate the Smart Code Key.

Alternative site testing-palm or forearm blood sampling

- To perform a test using samples obtained from alternative sites, install the clear cap on the lancing device (For more information on how to install, see the Instructions for the lancing device).
- To increase the blood flow, massage the puncture area of palm or forearm for a few seconds.
- Immediately after massaging the puncture area, press and hold the lancing device with the clear cap against palm or forearm.
- Then press the release button.
- Continue holding the lancing device against palm or forearm and gradually increase pressure for a few seconds until the blood sample size is sufficient (Refer to Instructions for the lancing device)



- 6) Touch and hold the drop to the edge of sample entry until you hear a " beep " and the view window is totally filled with blood. If the view window is not totally filled with blood or the test does not start. Then discard the test strip and repeat the test with a new test strip.
- 7) You will see the countdown mode on the screen. After 8 seconds, the test result appears. For more information about your test results, please see your User Manual.

- 8) Pull off the depth adjustable cap. Without touching the used disposable lancet, stick the lancet tip into the protective cover.
- 9) Hold the release button in one hand and pull on the plunger in the other hand will safely eject the used disposable lancet.
- 10) Discard the used disposable lancet into an appropriate puncture-proof or biohazard container.
- 11) Replace the depth adjustable cap after finishing the test.

For more information on how to use your meter and understand your test results, see the User Manual.

Test Result

- Blood glucose test results are shown on the meter as mg/dL or mmol/L, depending on which unit of measurement you have chosen. Consult your doctor before making any changes to your diabetes medication program.
- If your blood glucose result is unusually high or low, or if you question your results, repeat the test with a new test strip. You can also run a quality control test with your **Rightest™** Check Key and **Rightest™** Control Solutions to check your meter and strip. If the test result still remains unusually high or low, contact your doctor immediately.
- If you are experiencing symptoms that are not consistent with your blood glucose test results and you have followed all the instructions in this manual, contact your doctor immediately.
- The **Rightest™** Meter displays results between 10 and 600 mg/dL or 0.6 and 33.3 mmol/L. If your test result is below 10 mg/dL (0.6 mmol/L), " Lo " will appear on the screen. Please repeat your test again with another strip. If you still get a " Lo " result, you should immediately contact your doctor.
- If your test result is above 600 mg/dL (33.3 mmol/L), " Hi " will appear on the screen. Please repeat your test again with another strip. If you still get a " Hi " result, you should immediately contact your doctor.

Expected values⁽¹⁾

Fasting Blood Glucose	
GLUCOSE LEVEL	INDICATION
From 70 to 99 mg/dL (3.9 to 5.5 mmol/L)	Normal fasting glucose
From 100 to 125 mg/dL (5.6 to 6.9 mmol/L)	Pre-diabetes (Impaired fasting glucose)
126 mg/dL (7.0 mmol/L) and above on more than one testing occasion	Diabetes

Precautions


- Check the expiration date printed on the package every time you use the strip. Do not use expired test strips.
- Close the vial cap immediately after taking test strip out from the vial.
- Do not perform quality control test with expired control solution.
- Do not bend or twist the test strip. Damage of test strip may cause wrong result.
- Do not reuse test strips.
- Do not reuse lancets. Discard used lancets properly.
- Wait at least 30 minutes to perform a test if you have moved the meter to an area of different temperature.
- If you want to purchase a new control solution, please contact your authorized Bionime representative.

Warning

- Keep the test strips or vial cap away from children. They may cause a choking hazard. If a test strip or vial cap is swallowed, contact your doctor immediately.

Limitations

- Grossly lipemic (fatty) samples may interfere with some methodologies. To be aware of such interferences, patients under the supervision of their doctor should have baseline glucose values established by a clinical laboratory method prior to starting glucose monitoring at home. These baseline values should be checked periodically thereafter.

-  **NOTE**

 - Suggest not to use this meter close to source of strong electromagnetic radiation, to avoid interference with proper operation.
 - Suggest to keep meter free of dust, water or any liquid.

- Store the strips in the original capped vial at temperatures between 4 °C to 30 °C (39 to 86 °F) and relative humidity below 90 %. Do not freeze.
- Replace the vial cap immediately and close tightly after taking test strip out from the vial. Do not leave the cap of vial opened. If the strip is exposed in the air too long, it will absorb the moisture and cause wrong test result.
- Use test strips within 3 months after first opening.

The measurement range of the **Rightest™** System is 10 to 600 mg/dL or 0.6 to 33.3 mmol/L

Please refer to the Quality Control section of the User Manual.

For more information on error messages and trouble shooting, please refer to the Error Messages and Trouble Shooting section of the **Rightest™** User manual.

You may also contact customer service by calling local distributor or email to rightest@bionime.com.
(At all other times, you could contact your doctor for assistance)

Detection Principle ⁽²⁾

Performance Characteristics

Precision

The precision was evaluated including (i) venous whole blood sample (ii) 3 levels glucose control solution in period of 10 days, by 10 meters and 3 batches of strips.

(i) Venous whole blood sample:

Glucose levels	P-01	P-02	P-03	P-04	P-05
(1) Total test numbers (n)	300	300	300	300	300
(2) Mean mg/dL (mmol/L)	44.1 (2.4)	105.3 (5.9)	126.9 (7.0)	198.0 (11.0)	364.3 (20.2)
(3) SD mg/dL (mmol/L)	1.3 (0.07)	2.3 (0.13)	2.9 (0.16)	3.6 (0.20)	5.8 (0.32)
(4) CV (%)	2.9 %	2.2 %	2.3 %	1.8 %	1.6 %

(ii) Control solution:

Glucose levels	CS-L	CS-N	CS-H
(1) Total test numbers (n)	300	300	300
(2) Mean mg/dL (mmol/L)	64.0 (3.6)	107.7 (6.0)	293.0 (16.3)
(3) SD mg/dL (mmol/L)	1.8 (0.10)	2.7 (0.15)	4.0 (0.22)
(4) CV (%)	2.9 %	2.5 %	1.4 %

The accuracy of the **Rightest™** Blood Glucose Monitoring System was demonstrated by comparing whole blood (plasma equivalent) glucose values on the **Rightest™** meter with plasma glucose values on a lab instrument. A total of 106 patients were enrolled. Each patient collected and tested their own blood samples (from the fingertip, palm and forearm) using the **Rightest™** System. Then the blood samples were centrifuged immediately after collection to obtain plasma. Analyze the plasma by the lab instrument - YSI 2300. 99.7 % of **Rightest™** meter values were within ± 15 mg/dL (0.83 mmol/L) of the YSI values at glucose concentrations < 100 mg/dL (5.55 mmol/L) and within ± 15 % at glucose concentrations ≥ 100 mg/dL (5.55 mmol/L). The results and differences between the two methods, **Rightest™** System and YSI 2300 (as the reference method) are proved in the tables below.

Table 1: represents samples for glucose concentrations < 100 mg/dL (5.55 mmol/L).

Difference range in values between the YSI value and the Rightest™ meter value	The percent (and number) of samples of alternative site were the difference between the Rightest™ and the YSI value within the range shown in the side row.		
	Fingertip	Palm	Forearm
Within ± 5 mg/dL (0.28 mmol/L)	64.5 % (40/62)	71 % (44/62)	66.1 % (41/62)
Within ± 10 mg/dL (0.56 mmol/L)	91.9 % (57/62)	96.8 % (60/62)	100 % (62/62)
Within ± 15 mg/dL (0.83 mmol/L)	100 % (62/62)	100 % (62/62)	100 % (62/62)

Table 2: represents samples for glucose concentrations ≥ 100 mg/dL (5.55 mmol/L).

Bias range in values between the YSI value and the Rightest ™ meter value	The percent (and number) of samples of alternative site were the bias between the Rightest ™ and the YSI value within the range shown in the side row.			
	Fingertip	Palm	Forearm	
	Within ± 5%	46.7 % (70/150)	54 % (81/150)	51.3 % (77/150)
	Within ± 10%	82 % (123/150)	89.3 % (134/150)	83.3 % (125/150)
	Within ± 15%	100 % (150/150)	100 % (150/150)	98.7 % (148/150)








* Acceptance criteria in ISO 15197 : 2013 are that 95 % of all differences in glucose values should be within ± 15 mg/dL (0.83 mmol/L) at glucose concentrations < 100 mg/dL (5.55 mmol/L), and within ± 15 % at glucose concentrations ≥ 100 mg/dL (5.55 mmol/L).
Note: For glucose concentrations < 100 mg/dL (5.55 mmol/L), difference values are expressed in mg/dL (mmol/L), and for glucose concentrations ≥ 100 mg/dL (5.55 mmol/L), difference values are compared in percentage.

The following compounds may interfere with the glucose measurement at the concentrations listed:

Uric acid ≥ 10 mg/dL (0.59 mmol/L)	L-Dopa ≥ 3 mg/dL (0.15 mmol/L)
Ascorbic acid ≥ 5 mg/dL (0.28 mmol/L)	Dopamine HCl ≥ 2 mg/dL (0.11 mmol/L)
Glutathione reduced ≥ 70 mg/dL (2.28 mmol/L)	Hemoglobin $\geq 6,000$ mg/dL (0.94 mmol/L)

Each Blood Glucose Test Strip contains the following reagents:	Glucose Oxidase (GOD)	12.3 %	
Potassium Ferricyanide	43.3 %	Non-reactive Ingredients	44.4 %

- 1) Diabetes Information - American Association for Clinical Chemistry (AACC) [Electronic Version] Retrieved Dec 21, 2015 from www.labtestsonline.org/understanding/analytes/glucose/test.html
- 2) In Vitro Diagnostics in Diabetes : Meeting the Challenge. Clinical Chemistry 45:9, 1596-1601 (1999).

IVD	For in vitro diagnostic use	EC REP	EU Representative	LOT	Lot number
	CE Mark with number of Notified Body		For single use only		Manufacturer
	Consult instructions for use		Biological risks		Use by
 Store between temperature 4°C and 30°C (39°F and 86°F)					

Version: October 2016